
MAXIPACK

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ES PCD - PROCESS CHALLENGE DEVICE LOAD RELEASE TYPE 5 CHEMICAL INTEGRATOR

► PRODUCT DESCRIPTION

The **MAXIPACK CHALLENGE TEST LOAD RELEASE** was developed as an effective device to evaluate the performance of steam sterilization cycles and assist those responsible for MSC (Materials and Sterilization Center) in quality management in the reprocessing of medical materials and devices. The system composed of a barrier that offers resistance to the action of the sterilizing agent associated with 01 **MAXITEST Type 5** Chemical Integrator that evaluates all the critical parameters of the cycle (steam quality, temperature, and time).

Classified as a challenge test because it represents the most difficult load configuration to be sterilized. Indicated to monitor sterilization cycles from 121°C to 134°C in loads that do not contain implantable health products (ANMSA RDC 15:2012).

The **MAXIPACK CHALLENGE TEST LOAD RELEASE** is available in boxes of 12 units according to the AAMI/ISO 11140-1, AAMI/ISO 17665-1 and AAMI ST 79 Standards.

► COMPOSITION

The package is produced in impermeable paper with an internal set of porous sheets that act as an important barrier system to air removal and vapor penetration. Complies with the standard established by AAMI/ST79. It also has an external chemical process indicator Type 1 in color, 01 **MAXITEST Chemical Integrator - Type 5**, housed in the center of the package, made of paper strip (medical grade) laminated front and back, printed with a thermochromic chemical formula to steam. The product is called an integrator because it performs a complete and integrated assessment of the critical parameters of the sterilization cycle (steam quality, time, and temperature). For internal use and developed in accordance with the AAMI/ISO 11140-1 standard, it guarantees efficiency and safety in the application. Dimensions: 19mm (H) x 78mm (W).

The compact size (9,5 cm x 9,5 cm) of the **MAXIPACK CHALLENGE TEST LOAD RELEASE** makes storage easy and is made entirely of recyclable materials, endorsing MAXIMED's commitment to conscientious use of resources and preservation of the environment.

► TECHNICAL INFORMATION

According to the determinations contained in the AAMI/ISO 11140-1 standard, the Type 5 chemical integrators must have a performance equivalent to or superior to those presented by the biological indicators and also relate to 03 sterilization cycles (temperature x exposure time) between 121°C and 134°C. The **MAXITEST Chemical Integrator - Type 5** fully meets these performance requirements and includes the following sterilization cycles:

121°C	127°C	134°C
20 minutes	12 minutes	3,5 minutes

► INSTRUCTIONS FOR USE

1. Take a **MAXIPACK CHALLENGE TEST LOAD RELEASE** package and note the identification data: Auto-Reader Incubators number, cycle number and processing date.

2. Place the package next to the materials to be sterilized in an area considered critical (of greatest challenge) of the Auto-Reader Incubators, usually located near the door and next to the drain. Make sure it is in a horizontal position and with the face with the pink Type 1 chemical indicator facing up.

3. Start the sterilization cycle corresponding to the load to be processed. After the end of the cycle, remove the **MAXIPACK CHALLENGE TEST LOAD RELEASE** package from the Auto-Reader Incubators, respecting the cooling time recommended by the manufacturer.

4. Confirm that the Type 1 External Chemical Indicator contained in the package has changed from its pink (original) color to brown/black, confirming that it has been exposed to steam at 121°C or above.

Important: color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

5. Open the **MAXIPACK CHALLENGE TEST LOAD RELEASE** package through the perforation located on the back and remove the **MAXITEST Chemical Integrator - TYPE 5** from the inside of the package.

6. Verify that the TYPE 5 Chemical Integrator has entirely changed from pink to brown or black. Compare the change to the Interpretation Guide printed on the strip and this Manual.

IMPORTANT: ACCEPT THE CHEMICAL INTEGRATOR WITH RESULTS EQUAL TO OR DARKER THAN THE REFERENCE (IN TONES FROM BROWN TO BLACK COLOR).

NOTE: Black color means the process was successful. However, the exposure time exceeded that required for efficient sterilization. These thermochromic ink shade variations occur due to the heating period of the Auto-Reader Incubators (heating ramp) until reaching the sterilization temperature. This heating time exposes the chemical integrators to sterilizing temperatures for a time much longer than the time-calibrated in a resistometer (ISO 11140). Thus, variations in shades of brown and black are acceptable due to the non-uniformity of sterilization cycles and different load conditions.

7. If there is a failure in the process, the material must be processed again using a new **MAXIPACK CHALLENGE TEST LOAD RELEASE**.

8. Record the results according to current legislation and the institution's protocol.

CAUTION: Do not use the **MAXIPACK CHALLENGE TEST LOAD RELEASE** to monitor ethylene oxide sterilization cycles, dry heat, chemical vapors, radiation or other sterilization processes that do not use saturated steam as a sterilizing agent.

► VALIDITY/ANVISA REGISTRATION

The **MAXIPACK CHALLENGE TEST LOAD RELEASE** package is valid for 3 (three) years from the date of manufacture. Information about the lot, manufacturing date, and expiration date are printed on the packaging and the **MAXITEST Type 5** Chemical Integrator. Product exempt from registration with ANVISA, according to RDC 185/2001.

► STORAGE

Store away from light at a temperature between 15°C and 30°C and relative humidity between 35% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemical products.

► DISPOSAL

MAXIMED uses lead-free technologies and other environmentally harmful heavy materials and prides itself on responsibly sustainable practices in its facilities, processes, products, and services. These actions are part of the company's values, which are concerned with people's health, quality of life, and the future of our planet. In this way, all the components of the **MAXIPACK CHALLENGE TEST LOAD RELEASE** - as well as the complete line of MAXIMED products - are made with 100% recyclable materials, following the institutional plan of ecological design. After use, dispose of them following the sanitary regulations in your region.