
CLICKTEST^{20M}

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EN BIOLOGICAL INDICATORS - STEAM 20M (IBVF20M)

PRODUCT DESCRIPTION

The **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** was developed to monitor sterilization cycles from 121°C to 134°C in vacuum pump-assisted Auto-Reader Incubators, with air removal by gravity and immediate use cycle. It makes it possible to efficiently and safely prove the microbiological lethality of sterilization cycles.

Its innovative design dispenses glass components and, combined with an exclusive activation system, simplifies use and provides safe handling by preventing any risk of accidents.

The **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** is available in boxes of 50 units according to the AAMIISO 11138-1, AAMIISO 11138-3, and AAMIISO 11140-1 standards.

COMPOSITION

Each **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** is made with a thermoplastic tube, stainless steel ACS, special filter, adhesive label with process chemical indicator (Type 1), the minimum population of 1×10^6 of *Geobacillus stearothermophilus* spores (ATCC 7953) inoculated in paper strip (spore carrier). It also has a purple-colored growth medium contained in a casing and thermoplastic lid with holes, plus a vapor permeable barrier.

INSTRUCTIONS FOR USE

1. Identify the **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** by writing the sterilizer number, cycle/load number, and processing date on its label.
2. Pack the biological indicator in a standardized package according to recommended sterilization practices and institution validated.
3. Place the package next to the materials to be sterilized, in an area considered critical (of greatest challenge) of the Auto-Reader Incubators, usually located near the door and next to the drain.
4. Select and activate the sterilization cycle corresponding to the load to be processed. At the end of the cycle, wait for the cooling time recommended by the Auto-Reader Incubators manufacturer, open the chamber and carefully remove the biological indicator from the internal package.
5. Verify that the Type 1 external chemical indicator (contained on the label) has changed from its pink (original) color to brown or black, confirming that it has been exposed to vapor at 121°C or above.

Important: The color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

6. Activate the **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)**. To do this, hold the ampoule and press the lid until you hear a "click", breaking the retaining seal of the culture medium. Shake it until all the liquid (culture medium) has completely descended from the lid and comes into contact with the spore-bearing strip (see Quick Manual in this document).

7. Take another **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** that has **NOT** been subjected to the sterilization process to be used as a positive control, in order to validate the incubation conditions (temperature and time), the viability of the spores and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing batch. Repeat step 1 (identification) and step 6 (activation).

8. Verify that the **fluorescence technology** incubator operates between 58°C and 62°C and that the programmed time is compatible with the two 20-minute biological indicators.

SUPER-FAST READING RESULT - FLUORIMETRIC

Incubate the O2 biological indicators (processed and positive control) for 20 minutes. After the end of the incubation, observe:

a) if the incubator reader did **NOT** detect a fluorimetric reaction in the **processed indicator**, it characterizes a **negative result** and a successful sterilization cycle due to the inactivation of the spore population.

b) detecting fluorescence by the incubator means a **positive result** for the survival of microorganisms and a failed sterilization cycle.

9. Release the load/load for use only if the **processed Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** presents a result "Negative - Absence of Microorganisms" and the **control indicator** results "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that the incubation was adequate.

Note: if the processed biological indicator shows a result "Positive - Presence of microorganisms" and, consequently, an irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance must be notified to adopt the appropriate measures.

10. Detach the biological indicator labels (processed and positive control), attach them to the designated book/file and record the results obtained.

11. Dispose of the biological indicator tubes/ampoules following the sanitary regulations in your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators; or in cycles of 134°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All materials that make up the **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** are 100% recyclable, following the MAXXIMED institutional plan of ecological design.

NOTE: if necessary, it is possible to confirm the result obtained by rapid reading, complementing it with colorimetric technology (reading by changing the color of the culture medium) as follows:

RESULTS BY VISUAL READING

- **Do not perform steps 9, 10, and 11 above.**

- Check that the incubator is operating between 58°C and 62°C and incubate the O2 biological indicators for 24 hours to obtain the result by visual reading through colorimetric methodology (change of color of the culture medium). Since the presence of microorganisms can be detected before 24 hours, reading every 04 hours is recommended. Results:

a) maintenance of the purple color of the culture medium indicates the absence of viable microorganisms. Therefore, a successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from purple to gray/greenish and later to yellow, indicates the presence/survival of microorganisms. Therefore, sterilization cycle failure and **unsatisfactory results**.

- Perform steps 9, 10, and 11 above.

IMPORTANT: The super-fast 20-minute fluorimetric reading and the optional 24-hour colorimetric reading are correlated with the 07 days incubation according to the FDA's reduced incubation time protocol (FIT), and its results have a correlation of $> 0 =$ to 97% of cases.

CAUTION: Do not use the **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** to monitor ethylene oxide sterilization cycles, dry heat, plasma/hydrogen peroxide vaporization, low-temperature steam and formaldehyde, chemical vapors, radiation, or other sterilization processes that do not use saturated steam as a sterilizing agent. Do not reuse biological indicators.

VALIDITY / ANVISA REGISTRATION

The **Clicktest Super-Fast Fluorimetric Biological Indicator - IBVF20M (STEAM)** has validity of 02 years from the date of manufacture. Information about the lot and due date are printed on the packaging and on the biological indicator tube. Product exempt from registration with Anvisa, according to RDC 185/2001.

STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemical products.

IMPORTANT: always transport, store, activate, incubate and discard in an upright position (lid up) to prevent leakage after product activation.