
CLICKTEST^{20M}

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PRODUCT DESCRIPTION

The **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** was developed to monitor plasma sterilization or hydrogen peroxide vaporization processes. It makes it possible to prove, in an extremely fast, efficient, and safe way, the microbiological lethality of the sterilization cycles of thermosensitive materials and devices.

Its innovative design dispenses glass components and, combined with an exclusive activation system, simplifies use and provides safe handling by preventing any risk of accidents.

The **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** is available in boxes of 50 units and conforms to AAMI/ISO 11138-1 and AAMI/ISO 11140-1 standards.

COMPOSITION

Each **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** is made with thermoplastic tube, stainless steel ACS, special filter, adhesive label with process chemical indicator (Type 1), minimum population of 1×10^6 of *Geobacillus stearothermophilus* spores (ATCC 7953) inoculated in strip of paper (spore carrier). It also has a purple-colored culture medium contained in a wrapper and a thermoplastic lid with holes, in addition to a barrier permeable to the sterilizing agent.

INSTRUCTIONS FOR USE

1. Identify the **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** by writing the sterilizer number, cycle/load number, and processing date on its label.

2. Pack the biological indicator in a standardized package according to recommended sterilization practices and institution validated.

3. Place the package together with the materials to be sterilized in the Auto-Reader Incubators, in an area considered critical (of greatest challenge) to the action of the sterilizing agent. Consult the manufacturer's manual for a better understanding because depending on the cycle used, this area/position may vary.

4. Select and activate the sterilization cycle corresponding to the load to be processed. At the end of the cycle, and with the necessary safety precautions, remove the biological indicator from the interior of the standard package.

5. Verify that the Type 1 external chemical indicator (contained on the biological indicator label) has changed from its red (original) color to yellow, only confirming that there has been exposure to the sterilizing agent.

Important: The color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

6. Activate the **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)**. To do this, hold the ampoule and press the lid down **until you hear a "click"**, breaking the retaining seal of the culture medium. Shake it **until all the liquid** (culture medium) **has completely descended** from the lid and comes into contact with the spore-bearing strip (see Quick Manual in this document).

7. Take another **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)**, which **has NOT been subjected to the sterilization process**, to be used as a positive control, in order to validate the incubation conditions (temperature and time), the viability of the spores and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing batch. Repeat step 1 (identification) and step 6 (activation).

8. Verify that the **fluorescence technology** incubator is operating between 58°C and 62°C and that the programmed time is compatible with the 20-minute biological indicators.

SUPER-FAST READING RESULT - FLUORIMETRIC METHODOLOGY

Incubate the O2 biological indicators (processed and positive control) for 20 minutes. After the end of the incubation observe:

a) if the incubator reader did **NOT** detect a fluorimetric reaction in the **processed indicator**, it characterizes a **negative result** and a successful sterilization cycle due to the inactivation of the spore population.

b) fluorescence detection by the incubator means a **positive result** due to the survival of microorganisms and a failed sterilization cycle.

9. Release the load/load for use only if the processed **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** presents a result "Negative - Absence of microorganisms" and the **control indicator** results "Positive - Presence of microorganisms", confirming that the cycle of sterilization was satisfactory and that the incubation was adequate.

Note: if the processed biological indicator shows a result "Positive - Presence of microorganisms" and, consequently, an irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance must be notified to adopt the appropriate measures.

10. Detach the biological indicator labels (processed and positive control), attach them to the designated book/record and record the results obtained for traceability or follow the institution's protocol for quality management.

11. Dispose of the biological indicator tubes/ampoules following the sanitary regulations in your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes Auto-Reader Incubators; or in cycles of 134°C for 04 minutes in Auto-Reader Incubators assisted by a vacuum pump.

All materials that make up the **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** are 100% recyclable, complying with the MAXIMED institutional plan of ecological design.

IMPORTANT: the final reading of the result of the **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** occurs in 20 minutes and is correlated with the incubation of 07 days, according to the FDA reduced incubation time (FIT) protocol, and its results correlate of > or = in 97% of cases.

ATTENTION: use the **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** to exclusively monitor the **PLASMA STERILIZATION/ HYDROGEN PEROXIDE VAPORIZATION** cycles. Do not reuse biological indicators.

VALIDITY / ANVISA REGISTRATION

The **Clicktest Super-Fast Fluorimetric Plasma Biological Indicator – IBPF20M (VH202)** is valid for 02 years from the date of manufacture. Lot and expiration date information is printed on the package and tube of the biological indicator. Product exempt from registration with Anvisa, according to RDC 185/2001.

STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Do not freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemical products.

IMPORTANT: always transport, store, activate, incubate and discard in an upright position (lid up) to prevent leakage after product activation.