
# CLICKTEST 48H

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### ETHYLENE OXIDE 48H (IBEC 48H)

#### ► PRODUCT DESCRIPTION

The **Clicktest Biological Indicator - IBEC48H (ETO)** was developed to monitor sterilization processes by ethylene oxide gas at temperatures between 30°C and 65°C. It makes it possible to efficiently and safely prove the microbiological lethality of sterilization cycles.

Its innovative design dispenses glass components and, combined with a unique activation system, simplifies use and provides safe handling while preventing any risk of accidents.

The **Clicktest Biological Indicator - IBEC48H (ETO)** is available in boxes with 50 units according to the AAMI/ISO 11138-1, AAMI/ISO 11138-3 and AAMI/ISO 11140-1 standards.

#### ► COMPOSITION

Each **CLICKTEST BIOLOGICAL INDICATOR - IBEC48H (ETO)** is made with a thermoplastic tube, ACS in stainless steel, a special filter, adhesive label with chemical process indicator (Type 1), the minimum population of  $1 \times 10^8$  of *Bacillus atrophaeus* spores (ATCC 9372) inoculated in a spore-carrying strip. It also has a culture medium of bluish green color contained in casing and thermoplastic cover with holes, plus a permeable barrier to the sterilizing agent.

#### ► INSTRUCTIONS FOR USE

1. Identify the **Clicktest Biological Indicator - IBEC48H (ETO)** by writing on its label the sterilizer number, cycle/load number, and processing date.

2. Pack the biological indicator in a standardized package according to the best sterilization practices and institution validated.

3. Place the test package next to the materials to be sterilized, in an area considered critical to the action of the sterilizing agent, usually located in the center of the sterilization chamber or away from the gas injector.

4. Select and trigger the sterilization cycle corresponding to the load to be processed.

5. At the end of the cycle, wait for the time required for the aeration of the materials (process for the disposal of gas residues from the inside of the products/packaging), take the test package and carefully remove the biological indicator from its interior.

6. Check that the external chemical indicator (Type 1), contained on the label, has uniformly changed from its original red color to variations between orange and yellow colors, confirming that the sterilizing agent has been exposed.

**Note:** variations between the orange and yellow colors mean that the process was successful, however, the exposure time or the concentration of the sterilizing agent exceeded the minimum required to change the chemical indicator.

Thus, variations in shades of orange and yellow colors are acceptable, due to the non-uniformity of sterilization cycles.

**Important:** color change is not evidence that the sterilization cycle was successful. If there is no color change it will be necessary to review the sterilization process.

7. Activate the **Clicktest Biological Indicator - IBEC48H (ETO)**. To do this, hold the ampoule and press the lid down **until you hear a click**, breaking the retaining seal of the growing medium. Shake it **until all the liquid** (culture medium) **has completely descends** from the lid and contacts the spore-carrying strip (see Quick Manual contained in this document).

8. Take another **Clicktest Biological Indicator - IBEC48H (ETO)**, which has **NOT undergone the sterilization process**, to be used as a positive control to validate the incubation conditions (temperature and time), the viability of the spores, and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing lot.

9. Repeat step 1 (identification) and step 7 (activation).

#### RESULT BY VISUAL READING

10. Check that the incubator is operating at 37°C (+/- 2°C) and incubate the O2 biological indicators for a period of 48 hours to obtain the result by visual reading using colorimetric methodology (color change of the culture medium). Since microorganisms can be detected before 24 hours, it is recommended to perform readings every 04 hours. Results:

**a)** maintenance of the bluish green color of the culture medium indicates absence of viable microorganisms and, therefore, successful sterilization process and **satisfactory result**.

**b)** change in the color of the culture medium, from bluish green to yellow, indicates presence/survival of microorganisms and, therefore, failure in the sterilization cycle and **unsatisfactory result**.

11. Release the labels of biological indicators for use only if the **processed Clicktest Biological Indicator - IBEC48H (ETO)** results in "Negative - Absence of Microorganisms" and the **control indicator** results in "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that incubation was adequate.

**Note:** if the biological indicator processed points to the result "Positive - Presence of microorganisms" and, consequently, irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance should be notified for the adoption of the appropriate measures.

12. Detach the labels of biological indicators (processed and positive control), attach them to the designated book/medical record and record the respective results obtained.

13. Dispose of the tubes/ampoules of the biological indicators according to the sanitary regulations of your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators or cycles of 134°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All materials that make up the **Clicktest Biological Indicator - IBEC48H (ETO)** are 100% recyclable, following the MAXXIMED institutional plan of ecological design.

**IMPORTANT:** The 48-hour colorimetric reading is correlated with the incubation of 07 days, according to the FDA's reduced incubation time protocol (FIT).

**WARNING:** Use the **Clicktest Biological Indicator - IBEC48H (ETO)** to **exclusively** monitor **ETHYLENE OXIDE** gas sterilization cycles. Do not reuse biological indicators.

#### ► VALIDITY/ANVISA REGISTRATION

**Clicktest Biological Indicator - IBEC48H (ETO)** has validity of 02 years from the date of manufacture. Information about the lot and due date are printed on the packaging and on the biological indicator tube. Product exempt from registration with Anvisa, according to RDC 185/2001.

#### ► STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Don't freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents and other chemicals.

**IMPORTANT:** always transport, store, activate, incubate and discard in an upright position (lid up) to prevent leakage after product activation.