

BOWIE & DICK

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ES AIR REMOVAL TEST AND STEAM PENETRATION READY PACKAGE TYPE - 4KG OR 7KG

PRODUCT DESCRIPTION

The **BOWIE & DICK MAXXIMED** Test is classified as Type 2 by AAMI/ISO 11140-1: Chemical Indicators for Use in Specific Tests. It has the purpose of evaluating the correct functioning of steam at the Auto-Reader Incubator assisted by a vacuum pump and accurately identifies faults, such as inadequate air removal, air leaks, poor steam penetration, presence of non-condensable gases, low-quality steam, temperature and/or incorrect exposure time.

For daily use, it reads the three critical parameters of the sterilization process - steam quality, temperature and exposure time - in specific cycles at 121°C for 15 minutes or 134°C for 3.5 minutes.

The **BOWIE & DICK MAXXIMED** Test is available in versions BD7KG according to the AAMI/ISO 11140 part 4 and BD4KG with AAMI/ISO part 5.

COMPOSITION

Each **BOWIE & DICK MAXXIMED** Test pack is made of impermeable material, which acts as a barrier to the sterilizing agent and it contains a pink chemical process indicator (Type 1). There are sheets of porous paper inside it, which ensure a high level of barrier: an early warning sheet and a test sheet in the center, both impregnated with a non-toxic thermochromic chemical indicator, without the presence of lead and heavy metals.

Its compact size (9.5 cm x 9.5 cm) facilitates storage and endorses MAXXIMED's commitment to conscientious use of resources and preservation of the environment.

INSTRUCTIONS FOR USE

1. Separate and note externally in a **BOWIE & DICK MAXXIMED** Test package the identification data: Auto-Reader Incubators number, cycle number and date.

2. After the preheating cycle and with the Auto-Reader Incubators empty, using perforated support, place the **BOWIE & DICK MAXXIMED** Test in the area considered most inaccessible (critical) to saturated steam, generally located near the door and between 10cm and 20cm above the drain.

3. Make sure the test pack is in a horizontal position, with the side containing the pink Type 1 chemical indicator facing up.

4. Activate the BOWIE & DICK Test cycle and at the end of the cycle, immediately remove the test pack from the interior of the Auto-Reader Incubators with precautions, as it will be hot.

5. Verify that the Type 1 External Chemical Indicator has changed from pink to brown or black shades. Open the package and read the results of the alert sheet and test sheets (see Interpretation Guide in this manual).

Attention: consider as an accepted/approved result the uniform change from pink to brown or black shades. This situation occurs as a function of the conditions generated in the conditioning phase (vacuum) and/or the time spent in the heating ramp phase inside the Auto-Reader Incubators chamber (time that takes from the end of the vacuum phase until reaching the temperature of sterilization). These phases expose the BOWIE & DICK Test package to conditions that exceed those established by ISO 11140 and are possible to achieve only in resistors (BIER Vessel).

6. Approved result = correct functioning of the Auto-Reader Incubators and equipment released to start the daily sterilization routine.

Failed result = test cycle failure. Take a new **BOWIE & DICK MAXXIMED** Test pack and run a new cycle following all the previous steps. If the new test fails, interdict the equipment and call Clinical Engineering/Maintenance to solve the technical problem.

7. Make the appropriate notes on the test sheet and record the data in the Institution/MSD sterilization cycle book or control. The alert sheet and the test sheet are legal documents and must be filed in compliance with current legislation.

Note: the **BOWIE & DICK MAXXIMED** test is for single use and must not be reused under any circumstances.

STORAGE

Store away from light in a clean, dry, and ventilated place. Do not leave direct contact with the floor or sterilizing agents, and keep away from heat sources or products that can give off odors. Observe the temperature range between 15°C and 30°C and relative humidity between 30% and 60% for storage and transport. Do not freeze.

DISPOSAL

MAXXIMED uses lead-free technologies and other heavy materials harmful to the environment and health and prides itself on maintaining sustainable practices in its facilities, processes, products and services.

These practices are part of the company's values, which are concerned with people's health, quality of life and the future of our planet.

All components of the BOWIE & DICK test - as well as the complete line of MAXXIMED's products - are made with 100% recyclable materials, following the institutional plan of ecological design.

After use, dispose of it following the sanitary regulations in your region.

VALIDITY

The **BOWIE & DICK MAXXIMED** test is valid for 03 years from the date of manufacture. Lot and expiration date information is printed on the package, early warning sheet and test sheet.

NORMAS

- ▶ BD4KG according to the AAMI/ISO 11140-1 / AAMI/ISO 11140-5.
- ▶ BD7KG according to the AAMI/ISO 11140-1 / AAMI/ISO 11140-4.
- ▶ Product of mandatory use in daily routines to evaluate the operation of steam autoclaves assisted by a vacuum pump, according to article 93 of RDC 1/2012 ANMISA.
- ▶ Exempt from ANMISA registration as it is not considered a health product, according to RDC 185/2001.