
CLICKTEST^{24H}

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EN BIOLOGICAL INDICATORS LOW TEMPERATURE STEAM AND FORMALDEHYDE 2H (IBFF2H)

PRODUCT DESCRIPTION

The **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** was developed to monitor sterilization processes by mixing low temperature and formaldehyde steam at temperatures between 50°C and 80°C. It makes it possible to efficiently and safely prove the microbiological lethality of sterilization cycles.

Its innovative design dispenses glass components and, combined with a unique activation system, simplifies use and provides safe handling while preventing any risk of accidents.

The Clicktest Biological Indicator - IBFC24H (LTSF) is available in boxes with 50 units according to the AAMI/ISO 11138-1, AAMI/ISO 11138-5 and AAMI/ISO 11140-1 standards.

COMPOSITION

Each **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** is made with a thermoplastic tube, ACS in stainless steel, a special filter, adhesive label with chemical process indicator (Type 1), the minimum population of 1×10^9 of *Geobacillus stearothermophilus* spores (ATCC 7953) inoculated in a spore-carrying strip. It also has a purple culture medium contained in casing and thermoplastic cover with holes, plus a permeable barrier to the sterilizing agent.

INSTRUCTIONS FOR USE

1. Identify the **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** by writing the sterilizer number, cycle/load number, and processing date on your label.

2. Pack the biological indicator in a standardized package according to the best sterilization practices and institution validated.

3. Place the package next to the materials to be sterilized in an area considered critical (most challenging) of the Auto-Reader Incubators, usually located in the load center or near the door and next to the drain.

4. Select and trigger the sterilization cycle corresponding to the load to be processed. At the end of the cycle, following the safety recommendations of the Auto-Reader Incubators' manufacturer, open the chamber and carefully remove the biological indicator from the inside of the package.

5. Check that the external chemical indicator Type 1 (contained in the biological indicator label) has changed its red (original) color to green, confirming exposure to low temperature and formaldehyde steam.

Important: color change is not evidence that the sterilization cycle was successful. If there is no color change, it will be necessary to review the sterilization process.

6. Activate the **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)**. To do this, hold the ampoule and press the lid down **until you hear a click**, breaking the retaining seal of the growing medium. Shake it **until all the liquid** (culture medium) **has completely descended** from the lid and contacts the spore-carrying strip (take a look at the quick manual contained in this document).

7. Take another **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)**, which **has NOT been subjected to the sterilization process** to be used as a positive control to validate the incubation conditions (temperature and time), assess the viability of the spores and the ability of the culture medium to promote spore germination. Both biological indicators (processed and positive control) must belong to the same manufacturing lot. Repeat step 1 (identification) and step 6 (activation).

8. Check that the incubator operates between 58°C and 62°C and that the programmed time is compatible with the 24h biological indicators. Incubate the 02 biological indicators for 24 hours to obtain the result by visual reading using

colorimetric methodology (change of color of the culture medium). Considering that the presence of microorganisms can be detected before 24 hours, it is recommended to perform readings every 04 hours. Results:

a) maintenance of the purple color of the culture medium indicates the absence of viable microorganisms and, therefore, a successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from purple to gray/greenish and later to yellow, indicates the presence/survival of microorganisms and, therefore, failure in the sterilization cycle and **unsatisfactory result**.

9. Release the loading/load for use only if the processed **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** results in "Negative - Absence of Microorganisms" and the results in "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that incubation was adequate.

Note: if the biological indicator processed points to the result "Positive - Presence of microorganisms" and, consequently, irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance should be notified for the adoption of the appropriate measures.

10. Detach the labels of biological indicators (processed and positive control), attach them to the designated book/medical record and record the respective results obtained.

11. Dispose of the tubes/ampoules of the biological indicators according to the sanitary regulations of your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators or in cycles of 134°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All materials that make up the **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** are 100% recyclable, following the MAXIMED institutional plan of ecological design.

IMPORTANT: The 24-hour colorimetric reading is correlated with incubation of 07 (seven) days, according to the FDA's reduced incubation time protocol (FIT), and its results have a correlation of $> =$ to 97% of cases.

WARNING: Do not use the **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** to monitor sterilization cycles by ethylene oxide, dry heat, hydrogen peroxide, plasma/vaporization, saturated vapor, chemical vapors, radiation, or other sterilization processes that do not use low-temperature vapor and formaldehyde as a sterilizing agent. Do not reuse biological indicators.

VALIDITY/ANVISA REGISTRATION

The **Clicktest Colorimetric Biological Indicator - IBFC24H (LTSF)** has validity of 2 (two) years from the date of manufacture. Information about the lot and due date are printed on the packaging and on the biological indicator tube. Product exempt from registration with Anvisa, according to RDC 185/2001.

STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Don't freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents and other chemicals.

IMPORTANT: Always transport, store, activate, incubate and discard in the upright position (cover up) to prevent leaks after product activation.