
CLICKTEST 4H

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EN BIOLOGICAL INDICATORS ETHYLENE OXIDE 4H (IBEF 4H)

PRODUCT DESCRIPTION

The **Clicktest Biological Indicator - IBEF4H (ETO)** was developed to monitor sterilization processes by ethylene oxide gas at temperatures between 30°C and 65°C. It makes it possible to efficiently and safely prove the microbiological lethality of sterilization cycles.

Its innovative design dispenses glass components and, combined with a unique activation system simplifies use and provides safe handling while preventing any risk of accidents.

The **Clicktest Biological Indicator - IBEF4H (ETO)** is available in boxes with 50 units according to the AAMISO 11138-1, AAMISO 11138-3 and AAMISO 11140-1 standards.

COMPOSITION

Each **Clicktest Biological Indicator - IBEF4H (ETO)** is made with a thermoplastic tube, ACS in stainless steel, a special filter, adhesive label with chemical process indicator (Type 1), the minimum population of 1×10^8 bacillus atrophaeus spores (ATCC 9372) inoculated in a strip carrying spores. It also has a culture medium of bluish green color contained in casing and thermoplastic cover with holes, plus a permeable barrier to the sterilizing agent.

INSTRUCTIONS FOR USE

1. Identify the **Clicktest Biological Indicator - IBEF4H (ETO)** by writing on your label the sterilizer number, cycle/load number, and processing date.

2. Pack the biological indicator in a standardized package according to the best sterilization practices and validated in the institution.

3. Place the test package next to the materials to be sterilized in an area considered critical to the action of the sterilizing agent, usually located in the center of the sterilization chamber or away from the gas injector.

4. Select and trigger the sterilization cycle corresponding to the load to be processed.

5. At the end of the cycle, wait for the time required for the aeration of the materials (process for the disposal of gas residues from the inside of the products/packaging), take the test package and carefully remove the biological indicator from its interior.

6. Check that the external chemical indicator (Type 1), contained on the label, has uniformly changed from its original red color to variations between orange and yellow colors, confirming that the sterilizing agent has been exposed.

Note: variations between the orange and yellow colors mean that the process was successful, however, the exposure time or the concentration of the sterilizing agent exceeded the minimum required to change the chemical indicator.

Thus, variations in shades of orange and yellow colors are acceptable, due to the non-uniformity of sterilization cycles.

Important: color change is not evidence that the sterilization cycle was successful. If there is no color change it will be necessary to review the sterilization process.

7. Activate the **Clicktest Biological Indicator - IBEF4H (ETO)**. To do this, hold the ampoule and press the lid down **until you hear a click**, breaking the retaining seal of the growing medium. Shake it **until all the liquid** (culture medium) **has completely descended** from the lid and contacts the spore-carrying strip (take a look at the quick manual contained in this document).

8. Take another **Clicktest Biological Indicator - IBEF4H (ETO)**, which **has NOT undergone the sterilization process**, to be used as a positive control in order to validate the incubation conditions (temperature and time), the viability of the spores and the qualification of the culture medium. Both biological indicators (processed and positive control) must belong to the same manufacturing lot. **Repeat step 1 (identification) and step 7 (activation)**.

9. Check that the incubator with fluorescence technology is operating at 37°C ($\pm 2^\circ\text{C}$) and that the programmed time is compatible with that of the 04 hour biological indicators.

QUICK READING RESULT - FLUORIMETRIC METHODOLOGY

Incubate the 02 biological indicators (processed and positive control) for a period of 04 hours. After the end of incubation observe the results:

a) if the incubator reader did **NOT** detect fluorimetric reaction in the **processed indicator**, it characterizes **negative result** and successful sterilization cycle by inactivation of the spore population.

b) detection of fluorescence by the incubator means a **positive result** by the survival of microorganisms and **a failed sterilization cycle**.

10. Release the loading/load for use only if the **Clicktest Biological Indicator IBEF4H (ETO)** **processed** present results "Negative - Absence of microorganisms" and the **control indicator** result "Positive - Presence of microorganisms", confirming that the sterilization cycle was satisfactory and that the incubation and fluorescence detection were adequate.

Note: if the biological indicator processed points to the result "Positive - Presence of microorganisms" and, consequently, irregularity in the sterilization cycle, the Hospital Infection Control Commission and Engineering/Maintenance should be notified for the adoption of the appropriate measures.

11. Detach the labels of biological indicators (processed and positive control), attach them to the designated book/medical record and record the respective results obtained.

12. Dispose of the tubes/ampoules of the biological indicators according to the sanitary regulations of your region. Positive indicators can be sterilized in cycles of 121°C for 15 minutes or 132°C for 10 minutes in gravitational Auto-Reader Incubators, or cycles of 134°C for 04 minutes in vacuum pump-assisted Auto-Reader Incubators.

All materials that make up the **Clicktest Biological Indicator - IBEF4H (ETO)**, are 100% recyclable, obeying the MAXXIMED institutional plan of ecological design.

NOTE: if necessary it is possible to ratify the result obtained by quick reading complementing with colorimetric methodology (color change reading of the culture medium) as follows:

RESULT BY VISUAL READING

- Do not perform the previous steps 10, 11 and 12.

- Check that the incubator is operating at 37°C ($\pm 2^\circ\text{C}$). Incubate the 02 biological indicators for a period of 48 hours to obtain the result by visual reading through colorimetric methodology (change of color of the culture medium). Considering that the presence of microorganisms can be detected before 24 hours, it is recommended to read every 04 hours. Results:

a) maintenance of the bluish green color of the culture medium indicates absence of viable microorganisms and therefore, successful sterilization process and **satisfactory result**.

b) change in the color of the culture medium, from bluish green to yellow, indicates presence/survival of microorganisms, and therefore, failure in the sterilization cycle and **unsatisfactory result**.

- Perform the previous steps 10, 11 and 12.

IMPORTANT: The rapid fluorimetric reading of four 04 hours and the optional colorimetric reading of 48 hours are correlated with the incubation of 07 days, according to the FDA's reduced incubation time protocol (FIT).

WARNING: Use the **Clicktest Biological Indicator - IBEF4H (ETO)** to **exclusively** monitor **ETHYLENE OXIDE** gas sterilization cycles. Do not reuse biological indicators.

▶ VALIDITY/ ANVISA REGISTRATION

Clicktest Biological Indicator - IBEF4H (ETO) has a validity of 02 years from the date of manufacture. Information about the lot and due date is printed on the packaging and the biological indicator tube. Product exempt from registration with Anvisa, according to RDC 185/2001.

▶ STORAGE

Store away from light at temperatures between 15°C and 30°C and relative humidity between 35% and 60%. Don't freeze. Do not store in contact with floors and walls or near heat sources, sterilizing agents, and other chemicals.

IMPORTANT: always transport, store, activate, incubate and discard in an upright position (lid up) to prevent leakage after product activation.